

Cryogenic Sample Archiving Made Simple

Sample archiving package for cryogenic storage for only £ 1,875.

That's the same as saying "the laptop comes free!"



CONTAINS

Xtr-1elite Single Tube Reader

Cat. No. 20-1003. Single tube reader for all 2D data-matrix labelled tubes, hardware decoding with USB connection and keyboard wedge. List £1,250.

EasyTrack 2ml 2D DataMatrix-coded CryoTubes (5 packs) with reusable 48 format rack

Cat. No. 65-7501. 2 ml polypropylene, round bottom tubes with external thread, natural screw cap and starfoot. 0.4-1.8 ml working volume, 12.5 x 48 mm, natural, sterile, 100 per pack. List £97/pack.

Cat. No. 65-8501. 48 Way reusable rack for 2D-coded Cryotubes. Twist lock and open bottom for reading on Xtr range of rack reader. List £140.

Acer Netbook laptop

Xtr-1elite Single Tube Reader

Ideal for instant tube verification and initial ID entry.

The Xtr-1elite does not require any software installation or CPU usage and so can be used alongside old, or low specification, PCs.

Compatible with all versions of Windows™ operating systems.

Tube data can be entered into any application.

EasyTrack 2ml 2D dataMatrix-coded CryoTubes

Sterile polypropylene, round bottom tubes perfect for use in vapour phase liquid nitrogen storage conditions.

Unique 2D DataMatrix code permanently adhered to the base of each tube.

Automation friendly with external thread to reduce risk of contamination.

Starfoot design provides added stability. Free from DNase, RNase, DNA and endotoxins.

Acer Netbook. (Acer AOA110-Bb Atom N270 1GB 16GB XP Home Blue)

A powerful netbook in a diminutive form. Fast and simple, this incredibly mobile laptop is easily moved around sample storage facilities.

For further information, or to place an order, simply contact FluidX on:

Tel: +44 (0) 1625 861614

Email: order@fluidx.co.uk

Offer expires 30th June 2009.

Standard terms and conditions apply.

